# Patient Information for STAAR Surgical Implantable Collamer Lenses (ICLs)

## EN - English

### General Information

The STAAR Surgical Implantable Collamer Lenses (ICLs) are refractive lenses made from a material called Collamer®, a polymer proprietary to STAAR Surgical that is made of a similar material to that found in soft contact lenses combined with a small amount of purified porcine collagen. It also contains an ultraviolet light filter. Collamer is biocompatible (does not cause a reaction inside the eye) and is stable.

### Expected Device Lifetime

STAAR ICLs are long term implants. The long-term effect of the lens has not been determined.

### Patient Follow-up Information


Implantation of the ICL is a surgical procedure, and as such, carries potentially serious risks. The following represent potential complications/adverse events:

- under/over correction of vision,
- additional surgeries,
- ICL removal, replacement, or repositioning,
- cataract formation,
- temporary or permanent loss of best corrected vision,
- raised pressure inside the eye,
- a loss of cells on the innermost surface of the cornea,
- a reddening of the observable, white portion of the eyeball and inner eyelid,
- abnormal fluid build-up/swelling in the cornea,
- severe infection or inflammation of the entire eyeball,
- significant glare and/or halos around lights,
- blood in the eye,
- pus in the eye,
- eye infection,
- ICL dislocation,
- swelling in the area responsible for fine (reading) vision on the back surface of the eye,
- non-reactive pupil,
- the inability of fluid to flow from the back chamber of the eye to the front chamber frequently blocking drainage of fluid out of the eye and raising the pressure in the eye,
- severe inflammation of the eye, inflammation in the front chamber, middle layer of tissue, or other portion of the eye,
- loss of a clear gel-like material from the farthest back chamber of the eye during a surgical procedure,
- separation of the retina from its natural position on the back surface of the eyeball, and
- corneal transplant.

If you notice any sudden decrease in your vision, you should contact your doctor immediately. After ICL surgery it is important that you follow your doctor's recommendations for eye care and follow-up visits. The patient implant card should be used to provide information regarding your implanted ICL to any eye doctor that you visit.

### Patient Implant Card

The implant card provided to you by your surgeon contains information specific to the individual STAAR ICL that was implanted.

Information included on your implant card:

1. Your name, the date the ICL was implanted, and the name and address of the doctor
2. A description of the device that was implanted
3. The website to find important safety information about the ICL
4. The name, address, and website of the ICL manufacturer
5. The unique device identification (UDI) number and barcode. This number is specific to the manufacturer, manufacturing location, model, size, and power that was implanted
6. The date of surgery and name of the surgeon. These fields may be blank because the information is on the other side of the card.
7. The power, size of the optic and the overall size of the ICL that was implanted.
8. Model of ICL that was implanted and lens length
9. Serial number of the ICL that was implanted
10. Eye that was implanted. OD = right eye and OS = left eye.

### Manufacturer contact information

STAAR Surgical AG  
Hauptstrasse 104 CH-2560, Nidau  
Switzerland  
www.staar.com

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